



bioVis

VISCOELASTIC
SOLUTIONS
for intraocular use

MEDICENTUR



VISCOELASTIC SOLUTIONS FOR ALL KIND OF CATARACT SURGERIES

- ✓ EXCELLENT BIOCOMPATIBILITY
- ✓ PERFECT ENDOTHELIAL CELL PROTECTION
- ✓ OPTIMAL ANTERIOR CHAMBER AND CAPSULAR BAG MAINTENANCE
- ✓ EASY TO INJECT AND EASY TO REMOVE

PRODUCT	BioVis 1.6%	BioVis 1.8%	BioVis 3.0%
CONCENTRATION	16 mg/ml	18 mg/ml	30 mg/ml
STRUCTURE	Cohesive	Cohesive	Dispersive
ORIGIN	Biofermented	Biofermented	Biofermented
MOLECULAR WEIGHT	Appr. 3.0MDa	Appr. 3.0MDa	Appr. 0,75MDa
ZERO SHEAR VISCOSITY	400.000 mPas	600.000 mPas	30.000 mPas
OSMOLALITY	300 – 350 mOsm/kg	300 – 350 mOsm/kg	300 – 350 mOsm/kg
PH VALUE	6.8 – 7.6	6.8 – 7.6	6.8 – 7.6
VOLUME	1,1 ml	1,1 ml	1,1 ml
CANNULA	27G	27G	25G
STERILIZATION METHOD	Steam & ETO	Steam & ETO	Steam & ETO
STORAGE CONDITION	2-25° C	2-25° C	2-25° C



MC_BIOVIS_IS_EN v001



Mediontur
Medical Engineering Ltd
export@mediontur.com
www.mediontur.com

Head Office
Herceghalmi Road 1.
2072 Zsámbék
Hungary

International Office
Chemin des Aulx 18
1228 Plan-les-Ouates
Geneva / Switzerland

